

What's New

Changes in version 5.3.1



Office management system by Auditdata.
Version 5.3.1 released for Denmark, Norway, Sweden, and the UK.

What's New

Client

Client information

v 5.3.1

- Patient title is now a configurable lookup in Administration which you can edit and add to.

Custom questionnaires

Custom questionnaires

v 5.3.1

- Flexible Questionnaires can now be imported
- Flexible Questionnaires Slider works better with negative values

Letters, attachments and reports

User reports

v 5.3.1

- A new AQP report shows patients who have passed their 3 year review point

Measurement and clinical

NOAH

v 5.3.1

- Noah 4.5.1 including Noah Mobile. Noah Mobile supported for apps on mobile devices using HIMSA cloud. Customers must undertake data security review before using.
- Noah Audiogram Standard 502, will support Shadow Masking if the Noah module supports it.

Surgery – optional feature

v 5.3.1

- ABG calculation rounding is improved

Scheduling

Calendar

v 5.3.1

- Calendar now works with single click selection instead of double click in Dayview and Referral module

Referral pathway

v 5.3.1

- Paired outcomes can now be automated to appear three years or more after the initial outcome

Views and Navigation

Mouse scrolling and keyboard shortcuts

v 5.3.1

- Improved keyboard shortcuts

Navigation

v 5.3.1

- The ribbon state, whether maximized or minimized is saved after logout