

What's New

Changes in version 5.4.5



Office management system by Auditdata.

Version 5.4.5 released for Denmark, Norway, Sweden, Finland, Ireland and the UK.

What's New

Clinical records

Actions

- To save time when linking a Noah audiogram to an appointment manually only one linking dialog will be shown rather than the current two.

Journal

- We have starting using a new RTF control in the Journal which allows us to reinstate the Search feature, for both the display and edit windows. This control should also increase the population performance by 25%.

Letters, attachments and reports

Documents and attachments

- The 'Hearing Centre' variable for Word documents has now been outdated replaced by Location variables. The old variable will still work in existing reports.

Measurement and clinical

Audiometer and tympan interfaces

- We've made fixes to various bugs in tympanometry importing particularly for Titan.
- Tympanometry from Noah recorded in HIMSA format 501 can now be imported natively without losing the advanced format 501 features.

Auditory Implant Module – optional

- A new Implant surgery overview report allows export of initial surgery and switch on dates to Excel to help keep an overview of who should be billed for maintenance in any month and who should be transferring from paediatric to adult teams.
- All Implant preferences on the Pre-op Assessment form including Implant and Processor can now be transferred to the Client Information form.
- Custom Questionnaires can now be attached to the Pre-Op Assessment and Switch on screens, allowing you to use tests that you create yourself as the baseline and ongoing measurements.
- In the Custom Questionnaire designer you can select the fields that will display as summary results on the Pre-Op Assessment and Switch on forms.
- You can now open any test attached to the Pre-Op Assessment and Switch on screens directly from those screens.

-
- We've added various new tests and made changes to others commissioned by the SaLT and ToD members of the Bristol Impant teams;
 - We've added a comment field to the LING test.
 - We've added the EARS test.
 - We've added the eCaLL test.
 - We've added the GFTA and Khan-Lewis tests.
 - We've added the Symbolic Play test.
 - We've added a new option of how to record results on the CELF School Age test .
 - We've extended the milestone values in the NAMES 'Auditory Milestones - Nottingham' test.

NOAH

- The Noah version has been upgraded to the highest currently available version Noah 4.9.2.

Scheduling

Booking

- All restrictions and warnings on having users from different Depts in one appointment or time units from more than one Dept in one booking have been removed. This allows meetings of staff from different departments to be made more easily. Restrictions for different Locations are unchanged.
- We've made improvements to the speed of the save, update and delete transactions for Bookings to try to minimize the risk of locking errors.

System

Technical foundation

- We've made big changes to replace technical components from XSAL for tables, compression, tooltips, bitmap handling, zipping and encryption. These were no longer supported in the latest version of the Windows runtime that we use so replacing them allows us upgrade that in the near future.

System Administration

User setup, rights and maintenance

- To improve reliability and performance when creating Users we are changing the technology that user Roles are based on.
- User rights to tables, Document Types, Journal keywords and Extended rights will not be able to be customized, all rights changes will have to be made at the Role level.

-
- Existing Roles will be converted to the new Role type and Users will keep all of their existing rights.
 - Users who are not part of a Role or whose rights have been customized compared to a Role will be not be part of a Role after upgrade and cannot be modified until added to a Role.
 - A new Migration page in Administration will allow these users to be easily added to Roles.
 - At the same time as the technology change a big feature change is that users will be able to be added to multiple Roles simultaneously. Their rights will be an accumulation of the rights of their assigned Roles. As well as a useful feature in its' own right this provides compatibility with how Security Groups work in Windows Active Directory for future comptibility with that.