

What's New

Changes in version 5.3.0

AuditBase

Office management system by Auditdata.
Version 5.3.0 released for Denmark, Norway, Sweden, and the UK.

What's New

Client

Client merge

v 5.3.0

- Client Merge is now an Extended User Right in AuditBase rather than a System administrator function in Administration. The performance of the merge process has been improved to make it quicker and more reliable.

Postcode validation

v 5.3.0

- The postcode validation must now use the detailed validation which is set by default for all UK regions and can be constrained for the specific postcodes allowed in your region.

Clinical records

Journal

v 5.3.0

- Improved reliability of Journal module refreshing.

Communication

Task list – optional feature

v 5.3.0

- Users in the Task List are arranged by those from the currently selected Location first.

Hearing instrument related

Earmould dialogue

v 5.3.0

- The 'Sent to manufacturer' date added from the earmould module to the earmould list.
- A new Auto Journal Event added based on the completion of the 'Ready' box in the earmould module.

Stock

v 5.3.0

- The Stock Issue from Noah dialog can't be closed without adding an Auto Journal Event where users must explain the reason for cancelling the dialog without issuing the device.
- Serial numbered and quantity-controlled stock items can now be issued without indicating the ear side.
- The 'Warranty months' for devices can now support three digits, therefore warranty periods up to 999 months can be recorded.

Letters, attachments and reports

Letters and printing

v 5.3.0

- Details of both cancelled and a new appointment resulting from the cancellation available as Letter Variables.

User reports

v 5.3.0

- The default reports that are not available for your country, contain a 'Draft' watermark.

Measurement and clinical

AuditBase Viewer

v 5.3.0

- When AuditBase or the Viewer is opened from another application and the password has expired, the password change dialog will open immediately.

Graphical speech – optional feature

v 5.3.0

- The symbol mapping for importing Graphical Speech Test curves is much simplified.

Screening and Diagnostic ABR and OAE

v 5.3.0

- Insert or TDH phones recorded on the ABR/OAE Test form.
- ABR/OAE Diagnostic and screening results can be printed.

Surgery – optional feature

v 5.3.0

- An Audiogram can now be attached to more than one Surgery record, so it can be the post-op outcome for one Surgery and the pre-op starting point for the next.
- For Otosurgery and Implant surgery, the selection of Referring Agents now defaults to Surgeons. If a user is himself a surgeon, the Surgeon field can default to a Referring Agent associated with the current user.

Scheduling

Toolbar

v 5.3.0

- The 'Repair' button on the Noah error banner renamed to 'Reinitialize'.

Waiting list

v 5.3.0

- The breach date of the current Referral stage can be shown on the Waiting List.

System

Workstation and server environment

v 5.3.0

- Upgraded to Gupta runtime 6.2.4 for bug fixing.

System Administration

Controlling connections

v 5.3.0

- The last eight characters of the computer name are used for better unique identification of the PC in the Connection Monitor.

User setup, rights and maintenance

v 5.3.0

- Users can be associated with a Referring Agent.