

What's New

Changes in version 5.4.2



Office management system by Auditdata.

Version 5.4.2 released for Denmark, Norway, Sweden, Finland, Ireland and the UK.

What's New

Client

Client information

- Swedish and Norwegian temporary CPR numbers with 60/40 added to date supported

GDPR

- GDPR - right to be forgotten
 - Permanent deletion of identifying information in AuditBase
- GDPR - right to be informed
 - Reports can be generated showing all the information collected on a named patient.
- GDPR - right to data portability
 - Patient data can be exported in CSV format

Hearing instrument related

Stock

- Patient name can be shown in Stock Statistics for quantity controlled Accessories
- Show which staff member a device is lent to in Stock Statistics
- Fast accessory issue - quickly issue multiple articles in one step
 - Issue multiple accessories from different collections
 - The patient's most recently issued accessories
 - Accessories linked to current serial controlled devices or from the same manufacturer
 - Accessories from lists that you can build from those frequently issued or manually built
 - Using the traditional group and manufacturer selection
- Double click to add to the basket or use the number control and arrows
- Increase or decrease the quantities in the basket with + and -
- Instant stock quantity feedback and level control
- Single journal event
- Keyboard shortcut to open from anywhere
- Links to device issued from Noah
- Danish barcode support reestablished

Letters, attachments and reports

Documents and attachments

- Tympanograms are available as variables for Word Documents
- Default Norwegian Graphical Speech report is now printed on one page

Measurement and clinical

NOAH

- New Noah version 4.8 supporting format 500 for fitting modules
- Audiograms with historical dates can be saved from Noah measurement modules

Surgery – optional feature

- New 'Chorda tympani' combo box

Other features

Attention list

- Treatment Summaries submitted electronically can be filtered by date on Attention List

SE eFreeCard

- Support of Swedish online reimbursement solution (eFrikort)
- Unsuccessful submissions to online reimbursement solution (eFrikort) on Attention list
- Delete registrations to online reimbursement solution (eFrikort)
- Registrations to online reimbursement solution (eFrikort) included on posting

SE Self Check-in with visit fee

- Self Check-in for SE with visit fee, eCard, waiting room and directions support
- Support of Booking- / Self-check-in numbers
- Booking- /Self-check-in number included in caller and SMS notification
- Self-check-in access control (Self-check-in allowed / Not allowed)
- Prohibit self-check-in by patients with protected identity
- Prohibit self-check-in for specific symbols
- Client identity verification registration
- New Waiting room functions
- Support waiting room selection

-
- Default waiting room based on symbol and location
 - Visit location with name and address of waiting room available (e.g. for SMS notification)
 - Anonymous numbers for anonymous waiting room calls
 - Route description from present location
 - New printing possibilities (includes route description and anonymous visit number)
 - Additions to user interface only visible where Self-check-in is used
 - New tab with Self-check-in related functions on Booking dialogs
 - New columns with Self-check-in related information in Time plan and Day view dialogs
 - New section with Self-check-in related information in Search for Booking dialog
 - Additional Visit charge features
 - Register payments from Self-check-in system
 - Default Appointment fee (based on symbol) available on booking dialogs
 - Configurable set of payment methods (in Sweden the default set includes Swish)
 - Free of charge age ranges (e.g. younger than 19 and older than 85 are free of charge)
 - Check-in and payment location available in AuditBase (can be included on receipt)
 - Jump out to reception view of Self-check-in system without reentering patient

Scheduling

Booking

- New way to preview Resource availability
- New option to 'Suppress colours on the Resource screen' on site level and user level

Waiting list

- Waiting Lists can be exported to Excel

System Administration

User setup, rights and maintenance

- Search for a User in User Administration and Security