

# Auditbase

Office Management System

Version 6.7.0 released for Denmark, Norway,  
Sweden, Finland, Ireland and the UK

## What's New

Changes in version 6.7.0

---

**Auditdata**

## What's New

---

# Audiogram

## Audiogram

- We added new variables to the 'Large audiogram printout' report template, including speech data and masking levels.

# Client

## Client search

- On client search using Netherlands Employee ID, the system adds leading zeros allowing the user to search more quickly and easily.

# Clinical records

## Auto Journal Events

- The Automated Journal Event for audiogram creation has now been extended to include when the audiogram is created in Noah as well as in Auditbase.

## Journal

- Stored procedure support has been improved for accessing Journal records less than 32kb quickly.
- Since there have been incidents of Journal corruption caused by pasting content that includes images into the Journal and images cannot be viewed in the Journal, users will now be prohibited from doing this.

# Documentation

## Documentation

- The Auditbase 'About box' has been updated for MDR requirements, including a
-

[link to the Instructions For Use on the Auditdata.com website.](#)

# Hearing instrument related

## Stock

- You can now register articles in Article Groups in the Register for Stock dialog for serial controlled devices, speeding up the process of registering new articles.
- The Supplier Article Number has now been added to the Register for Stock dialog so that users can verify that they are registering the correct article during the process.
- The Requisition dialogs and processes have undergone a series of improvements.
- The size of the Requisitions dialog can now be expanded by the user and the size will be memorized for each user.
- The Requisitions and Purchase orders dialogs now show the full name of the patient for whom a device has been ordered.
- Devices ordered for a patient can now be viewed on a new 'Requested' tab on the Device list showing their order status.
- The content of the Requisitions list can now be exported to CSV by users with Extended Rights enabling it to be used more easily with external ordering systems.
- The default size of the Device List dialog has been increased. It can now be expanded by the user, and the size will be memorized for each user.
- You can now record the batch/lot number of serial-controlled devices, either through barcode scanning or manually.
- The batch/lot number of devices is shown in the device details, the device history, and the client device list.
- You can search for devices by their batch/lot number and find devices in stock or issued to patients. Where they are in stock, you can make bulk changes to their stock status. Where they have been issued to patients, you can add those to a Client List or Onward Referrals and Reporting List for potential recall.
- Purchase Orders listed can now be filtered by Location.

## Interfacing

### Electronic submission

- The audiogram legend in Electronic Submission audiograms is now configurable

to prevent duplication of the legend where it already exists in the report template.

# Letters, attachments and reports

## Documents and attachments

- Having macros in Word templates may cause unreliability when working with Word. Since some customers rely on it heavily, we have made it optional now, so if needed, it can be enabled at the risk of some adverse effects on stability.

## User reports

- A new version of the BCIG Cochlear Referral report includes an age parameter allowing you to report on adults and children separately or together.

# Logging

## Logging

- The Auditbase error log can now be exported directly from Administration, allowing IT departments to review it more easily and share it with Auditdata support.

# Measurement and clinical

## Audiogram

- The Audiogram means control has changed. Now the menu item only opens the floating mean display dialog. The choice of which mean is displayed in the Information tab is now made in that tab and returns to the default when changing the client.
- Numerical speech data from the Speech tab can now be included on the Compare Audiograms dialog with an option in Administration.
- The Compare audiograms printout has also been extended to include speech data where required.
- Numerical speech data from the Speech tab can now be included on the Audiogram preview dialog with an option in Administration.

# NOAH

---

- The latest Noah engine version 4.15 is now certified for use in Auditbase. The latest Noah Audiogram Module 1.5.0 is included in the installation.
- The Noah session list is now more user friendly as it lists the module name and when that module is installed it shows the module icon next to each action.
- Separate components of module actions are now linked by ActionGroupGUID allowing any action component to be used to open the relevant parent action.
- Usernames no longer have to be unique to overcome an issue with Noah. A new 'Noah unique Full name' is now generated based on the existing username, truncated to the maximum length that Noah allows, and made unique with a numerical postfix.

## Surgery

- New fields for Drain, Facial nerve and Lateral semicircular canal have been added to the Otosurgery Surgery screen.
- All fields for all tabs in both Otosurgery and Stapes Surgery dialogs now have value sets that are configurable in Administration.
- Surgery reports now have an additional parameter allowing users to include only patients from a chosen Client List or Onward Referral and Reporting List.
- A new report 'Cholesteatoma, graft technique & materials according to drum state' has been added.

## Other features

### Laptop synchronization

- Laptop synchronisation has now been updated to allow the second part of the COSI and Glasgow Questionnaires to be completed on synchronising laptops when the first part was done in the central database.

## Scheduling

### Booking

- The export of a user's Resource calendar to Outlook can be disabled in Administration, although when enabled, the export can no longer overwrite existing Outlook bookings and does not contain patient names or other details.
- The appointment DNA invoice event now records the Appointment ID to correct an issue with the automated Invoice Server.

# Security

## Access security

- Default password security settings for new databases have been increased.

## Vulnerability reduction

- A potential vulnerability in Custom Reports has been eradicated through a feature change.
- The TeamViewer QuickSupport executable has been removed from the installation process since it has been associated with previous security vulnerabilities. The link from the Help menu will instead take the user to the URL to download the very latest version for installation.

# System Administration

## User setup, rights and maintenance

- The Extended Right "Set/remove stock item Recycled" has been renamed to "Remove stock item Recycled".