

Auditbase

Office Management System

Version 6.3.0 released for Denmark, Norway,
Sweden, Finland, Ireland and the UK

What's New

Changes in version 6.3.0

Auditdata

What's New

Client

Client information

- We have added support for the Swedish National Reserve Identifier
- We have added new configuration options to define what can be viewed and what is mandatory, to be set at the site and individual user levels.
 - Hide unwanted tabs in Client Information module
 - Hide unwanted fields in the Client Information Details tab
- Mandatory Client fields - as well as last name which is already mandatory now first name and date-of-birth can be required
- We've added a new tab for 'Record info' to show when a client's record was created and whether they originated from a merge or conversion
- New gender values of Non-Binary and Not applicable have been added to the list of available values
- We have a new customer in the Netherlands so have introduced the Dutch national ID setup and partial translation to Dutch language

Clinical records

Tinnitus Functional Index

- We have added a 'Stage' dropdown in the Tinnitus Functional Index using the same values as in the Tinnitus Handicap Index so you can record pre and post treatment stages specifically and report on them

Custom questionnaires

Custom questionnaires

- "Age at test" which has always appeared on Custom Questionnaires can now be excluded from a Questionnaire design

Letters, attachments and reports

Letters and printing

- TLS1.2 security level is now supported for emailing allowing the latest Office 365 to be used as the email server

User reports

- We've introduced better separator handling for export reports to Crystal Format
- For administrators managing Ear surgeries we have introduced a new user report to make it easier to overview where patients are on their progress through assessment and treatment stages.

Measurement and clinical

Audiogram

- Audiogram static norm curves can now be defined in Auditbase Administration, a specific curve can be selected for a specific workstation
- Based on data from ISO 7029 you can now display age and gender specific dynamic hearing loss distribution percentile curves. The percentile value used can be chosen in Auditbase Administration on a site level
- We have added a choice of whether male, female or no percentiles should be shown in the Audiogram module when on a specific client. This allows a client with one stated gender to use percentile values calculated for their sex at birth.
- We have added new configuration options to define what can be viewed and what is mandatory, to be set at the site and individual user levels.
- Hide any of the tabs in Audiogram module, except 'Information'
- Set the 'Audiogram type' as a mandatory field for saving a record

Audiometer and tympanometry interfaces

- Tympanometry ID can now be edited manually to make it possible to record the Tympanometer details when using different tympanometers in the same room.
- Where a Noah Tympanometry module requires its curve compensation can now be applied even when it is declared as already applied in the measurement.

AuditBase Screener

- New Hearing Screener product using auto mode in audiometer panels with a specific selection of modules and a bespoke navigation

NOAH

- We have now made it possible to install Auditbase alongside Noah System without interacting with it for PCs where Auditbase is only wanted for administrative work and Noah System is used on a separate set of patients

Scheduling

Booking

- When exporting the Auditbase schedule to Outlook for a Resource we are now removing patient names for information security. To avoid losing Outlook data we have removed the option to delete existing Outlook appointments.